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SMART EYE MASSAGER
Instruction manual

Smart eye massager is an innovative device with heating and vibration functions that provides relaxation, reduces fatigue and pressure in temples. The Bluetooth function gives you the ability to connect to your phone and listen to your favorite music while relaxing. The device is easy to use, with a clear display indicating battery level, vibration or heating modes. The massager is recommended for people who spend long hours in front of a computer, as well as those who suffer from migraines or wear glasses. The modern design, folding function and compact size make the device easy to take on the go.

Rules for safe use:

1. Read the instructions carefully before using the product for the first time.
2. The device should not be used if there are any contraindications. If in doubt, consult your doctor.
3. Always start massaging with the device from the lowest intensity level.
4. If you experience pain or discomfort while using the massager, turn off the device immediately and consult a doctor.
5. The massager should not be used while driving a car.
6. The device should not be used while bathing or in direct contact with water - this could cause an electric shock.
7. Do not use the device if damaged, and contact the service center.
8. It is not recommended for children to use the massager.
9. The massager should not be used while sleeping.

Contraindications:

The device should not be used by people:

- with eye diseases,
- after eye surgeries, operations,
- with glaucoma, cataracts, retinal detachment,
- with wounds, inflammation,
- after strokes,
- with hypertension, hypotension,
- sensitive to heat.

If you have any doubts about the use of the device, we recommend consulting your doctor.

How to use the eye massager?

1. Press and hold the power button to start the device. The heating function and Bluetooth connection will start automatically.
2. Another press of the button will switch to a higher heating level.
3. Another two quick presses of the button will switch to vibration and pressure boost mode. One click changes the intensity.
4. Another double press of the button starts both vibration and heating functions at the same time. A single touch of the button increases the intensity of both functions.
5. Another double press returns to the initial settings.
6. Press and hold the button to turn off the device.

Troubleshooting:

The device does not start:	<ul style="list-style-type: none">• make sure the device is on and that the battery is charged - if not, charge the device.
Device shuts down automatically:	<ul style="list-style-type: none">• possible low battery level. Charge the device,• device still does not start - contact service.
Unable to connect to bluetooth:	<ul style="list-style-type: none">• Turn the device off and on again, try to connect to the phone.• Make sure the device is not connected to another phone.

Included:

- Smart eye massager
- charge cable
- user manual